

LymphoSure

Lot Number: Syn270208A Normal

Expiry Date: 31-03-2008

Expected Values*				
	% Positive	Range (%)	Absolute number (/uL)	Range (/uL)
CD3+	73	59-88	1008	806-1209
CD3+CD4+	45	36-54	618	494-741
CD3+CD8+	25	20-31	350	280-419
CD3-CD16+56+	12	10-14	165	132-198
CD19+	10	8-12	138	110-165
Lymphocytes:	34	27-41	N/A	N/A
Monocytes:	4	3-5	N/A	N/A
Granulocytes:	56	44-67	N/A	N/A

* These values have been determined using Becton Dickinson MultiTest™ reagents with TruCount™ tubes, following the lyse, no-wash method.