

Setting a new standard with LymphoSure



LymphoSure is a new standard for flow cytometric analysis of lymphocyte subsets. It was developed by Synexa Life Sciences to overcome the inconsistent scatter properties and poor shelf life encountered with existing commercial standards.

LymphoSure can be used as an IQA and for qualifying reagents in a flow cytometry lab. It is compatible with the full range of Becton Dickinson equipment (FACSCalibur, FACSCount and FACScan).

LymphoSure is provided in 2.5ml EDTA tubes, in both Normal and Low CD4 ranges. Larger sizes are available by request.

LymphoSure offers a number of key advantages:

- Stability** Guaranteed to be stable until its expiry date, provided it is stored correctly.
- Long life** Shelf life of 30 days in customer's hands.
- Quality** Clear and consistent scatter properties, enabling high quality gating.
- Price** LymphoSure is priced very competitively at approximately \$100 per 2.5 ml pack. Prices may vary slightly by region.
- Availability** LymphoSure is produced in frequent, small batches, ensuring stability and availability.

For further information on LymphoSure, contact Paul O'Riordan at poriordan@synexagroup.com or call +27 21 933 9580.